

Instruments

VISIT...

LANZAROTE
Caliente.COM

Scalpel





Non toothed forceps



Toothed forceps

Holding skin



Skin retractor





retractor
during appendectomy

Morris retractor

Retraction of abdominal wall

Deaver's retractor



Self retaining abdominal retractor



Curved artery forceps


used to hold tissue & catch bleeding point



Mosquito forceps

used to hold fine bleeding point






Allis forceps



Babcock forceps

used to hold intestine



Ring forceps



Needle holder



Handling forceps



Crushing intestinal clamp

applied over the resected part of intestine




Non –crushing intestinal clamp

used to prevent spilling of intestinal contents during intestinal anastomosis



Vascular (Vena cava) clamp

can be used for complete vessel clamping or to isolate a part of the vessel wall by side clamping



Intestinal holding

e.g Appendicular holding





Proctoscope ano-rectal examination



Curette




Self retaining urinary catheter



Disposable Urinary catheter

**Nelaton
catheter**





Urine collection bag



Supra-pubic catheter



Stone forceps


Stone kidney




**Stone urinary
bladder**



Stone ureter




Ureteric catheter




Urethral dilator




Fogarty catheter



Naso-gastric (Ryle) tube



The marks on the catheter guide to the length which has been introduced from the nose to reach the stomach



Air way

to maintain upper
airway patency &
prevents the tongue
from falling backwards
and obstructing the
airway



Tracheostomy tube




Closed system suction drain


It is more blessed to give than to receive

Suction unit



Corrugated rubber drain





“T” tube (for drainage of the common bile duct if it has been opened and explored)





I.V. cannula

The color of the cap indicates the size



Transfusion set






Sengstaken –Blakemore tube

Temponade & gastric wash

INDEX



Bone cutting



Orthopedic instruments

Osteotome & Chesil



Periosteal elevator




Rib raspatory

used to separate the periosteum during rib resection



It is more blessed to give than to receive.

Giggly saw



Verrous needle

inflation of abdomen during laparoscopy



Ports of laparoscope for introduction of instruments




Maryland forceps

Used for dissection



Laparoscopic scissor



THANK YOU

Dr. WAEEL METWALLY

